-online.com genomics

Datasheet for ABIN4843036 Human VCX3B cDNA Clone in Bacterial Expression Vector (His-GST)

Overview

Quantity:	500 ng
Gene:	VCX3B
Species:	Human
Fusion tag:	His-GST
Insert:	cDNA
Vector:	Bacterial Expression Vector
Application:	Cloning (Clon)

Product Details

Purpose:Bacterial expression of Human VCX3B with His-GSTInsert Length:645 bpVector Backbone:pPB-His-GSTPromoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-XholFusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Vector Backbone:pPB-His-GSTPromoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Promoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
	e-S-Transferase Protein
Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminato	r primer: 5'-

Target Details

Gene:

VCX3B

Target Details	
Alternative Name:	VCX3B
Application Details	
Application Notes:	 The pPB vectors are low-medium copy number vectors in which the gene expression is driven by the strong T7 promoter. Below are some basic guidelines for using the pPB vectors for protein production: The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the host E. coli cell has a source of T7 RNA polymerase. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl β-D-1-thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM. The ideal concentration of IPTG must be determined empirically for each recombinant protein/cell-line. Similarly, the length of time and temperature for induction provide other variables that need to be optimized on a case-to-case basis. For toxic proteins, it is recommended to go for shorter induction time and also to try and suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins. Once grown for the desired length of time, harvest cells by centrifugation and either freeze the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the purification.
Restrictions:	For Research Use only
Handling	
Format:	Liquid
Buffer:	10 mM Tris-HCl, 1 mM EDTA, pH 8.0
Storage:	-20 °C

Expiry Date:

12 months

Publications

Storage Comment:

Product cited in:

Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)

1 year when stored at -20° C or lower in a non-frost free freezer.