

Datasheet for ABIN4932305

Human C20ORF202 ORF Clone in Mammalian Expression Vector (DYKDDDDK Tag)

Overview

Quantity:	10 µg
Gene:	C20ORF202
Species:	Human
Fusion tag:	DYKDDDDK Tag
Insert:	ORF
Vector:	Mammalian Expression Vector
Application:	Protein Expression (PEXP)

Product Details

Purpose:	Expression/transfection ready cDNA ORF clone of Human C20orf202 with C terminal DYKDDDDK tag is ideal for express proteins in E.coli & mammalian cells.
Brand:	GenEZ™
Insert Length:	369 bp
Vector Backbone:	pcDNA3.1+C-(K)-DYK
Promoter:	CMV Promoter
Selectable Marker:	Neomycin
Bacterial Resistance:	Ampicillin
Expression Type:	Transient, Stable
Sequence:	ATGTACAAGT CAAAGATCCC TCGGGCCCAG AACCAGGTCA GTGTCAAGGT CACTCCTAAG AACACTGAGA TGAAAATAGC AGAAGAGCCC AGCCCGAGTC TTGGGCAGAC CTTGGAGTGG CTGAGAAAGG AGCTGTCTGA GATGCAGATT CAAGATCAGA GTCTCCTGCT CAACTGAGG CATCTTCACA GTGTCCTGGA GGAGCTGCGT GCGGACAGCG CCCACTGGGA GGACGCCAGG

Order at www.genomics-online.com

USA & Canada: +1 877 302 8632 | support@antibodies-online.com

Product Details

TCCAGCGGAG GGACATCCCC CATCAGAGCT CGAGCGGGCT CCGAAGGCAG GGGCTGCCAG
CCGGTTTGCT CAAGGGGTCT GGCCCAGCTC CTCCGAGGGG AAGACAGCAG ACGAAGCTCT
CTCCCTTAA

Specificity: ORF Insert Method: CloneEZ® Seamless cloning technology, recombination-based cloning technology

Characteristics: Gene cDNA ORF clone sequences were retrieved from the NCBI Reference Sequence Database (RefSeq). These sequences represent the protein coding region of the gene cDNA ORF which is encoded by the open reading frame (ORF) sequence.

Sequencing Primer:

- Forward primer: 5'-TAATACGACTCACTATAGGG-3'
- Reverse primer: 5'-CCTCGACTGTGCCTTCTA-3'

Grade: End-sequenced

Components: The GenEZ ORF clone is delivered as 10 µg of lyophilized plasmid DNA in a vial.

Target Details

Gene: C20ORF202

Alternative Name: C20orf202 ([C20ORF202 Products](#))

Gene ID: 400831

NCBI Accession: [NM_001009612](#)

Application Details

Restrictions: For Research Use only

Handling

Format: Lyophilized

Storage: RT/-20 °C

Storage Comment:

- Keep the vial sealed and store at -20°C for long-term storage.
- Before use, centrifuge the vial at 6,000 g x g for 1 minute at 4°C.
- Open the lid and add 100 µl (or other volume depending on your desired final concentration) of distilled water (or TE buffer) to dissolve the DNA.
- If necessary, heat the solution at 50°C for 15 minutes to dissolve the DNA.
- Close the lid and vortex the vial for 1 minute.
- Aliquot the dissolved plasmid DNA and store in small aliquots at -20°C.

Order at www.genomics-online.com

USA & Canada: +1 877 302 8632 | support@antibodies-online.com

Handling

Expiry Date: 12 months

Publications

Product cited in: Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)