

Datasheet for ABIN5519539

TransStart® Green qPCR SuperMix

Overview

Quantity:	1 mL
Application:	Quantitative real-time PCR (qPCR)

Product Details

Purpose:	TransStart® Green qPCR SuperMix is a ready-to-use qPCR cocktail containing all components, except primer and template.
Brand:	TransStart®
Specificity:	TransStart® Taq DNA Polymerase, hot start with double blocking technique, improves sensitivity, enhances specificity and generates more accurate data. Double cation (K ⁺ , NH ₄ ⁺) buffer enhances specificity and reduces primer-dimer formation. Passive reference dyes are provided for different qPCR instruments.
Characteristics:	It contains TransStart® Taq DNA Polymerase, SYBR Green I, dNTPs, PCR enhancer and stabilizer. qPCR SuperMix is provided at 2x concentration and can be used at 1x concentration by adding template, primer, passive reference dye (optional) and ddH ₂ O.

Application Details

Application Notes:	Optimal working dilution should be determined by the investigator.
Restrictions:	For Research Use only

Handling

Storage:	-20 °C
Storage Comment:	at -20°C for two years
Expiry Date:	24 months

Order at www.genomics-online.com

USA & Canada: +1 877 302 8632 | support@antibodies-online.com

Publications

Product cited in: Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)