-online.com

Datasheet for ABIN3234429 Human FAM90A14P CRISPR gRNA + Cas9 in Mammalian Expression Vector (Myc-DYKDDDDK Tag)

Overview

| Quantity: | 1 kit |
|-----------------------|--|
| Gene: | FAM90A14 (FAM90A14P) |
| Species: | Human |
| Fusion tag: | Myc-DYKDDDDK Tag |
| Insert: | gRNA + Cas9 |
| Vector: | Mammalian Expression Vector |
| Application: | Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Knockout Kit for Human FAM90A14 via CRISPR. |
| Vector Backbone: | pCas-Guide |
| Promoter: | U6 Promoter, Enhanced CMV Promoter |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Transient |
| Characteristics: | The FAM90A14 kit is designed based on the best knowledge of CRISPR technology. The system has been functionally validated for knocking-in the cassette downstream the native promoter. The efficiency of the knock-out varies due to the nature of the biology and the complexity of the experimental process. |
| Sequencing Primer: | CF3 (ACGATACAAGGCTGTTAGAGAG) |
| Components: | FAM90A14 gRNA vector 1 in pCAS-Guide vector. FAM90A14 gRNA vector 2 in pCAS-Guide vector. Donor vector containing Left and right homologous arms and GEP-Puro functional cassette. |

• Scramble sequence in pCas-Guide vector

Target Details

| Gene: | FAM90A14 (FAM90A14P) |
|---------------------|---|
| Alternative Name: | FAM90A14 |
| Application Details | |
| Application Notes: | Knock-in GFP reporter for promoter study.Knock-out genes at chromosomal level. |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Lyophilized |
| Storage: | -20 °C |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991) |