-online.com genomics



Datasheet for ABIN3392430

Human SPDYE2B cDNA Clone in Mammalian Expression Vector

protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Overview	
Species: Human Insert: cDNA Vector: Mammalian Expression Vector Application: Protein Expression (PExp) Product Details Purpose: Untagged full-length cDNA clone from Human SPDYE2L is ideal for over-expression of na protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Quantity:	10 μg
Insert: cDNA Vector: Mammalian Expression Vector Application: Protein Expression (PExp) Product Details Purpose: Untagged full-length cDNA clone from Human SPDYE2L is ideal for over-expression of na protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Gene:	SPDYE2B
Vector: Application: Protein Expression (PExp) Product Details Purpose: Untagged full-length cDNA clone from Human SPDYE2L is ideal for over-expression of na protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. These cDNA clones are ideal for over-expression of native proteins for functional studies.	Species:	Human
Application: Protein Expression (PExp) Product Details Purpose: Untagged full-length cDNA clone from Human SPDYE2L is ideal for over-expression of na protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	nsert:	cDNA
Product Details Purpose: Untagged full-length cDNA clone from Human SPDYE2L is ideal for over-expression of na protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Vector:	Mammalian Expression Vector
Purpose: Untagged full-length cDNA clone from Human SPDYE2L is ideal for over-expression of na protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Application:	Protein Expression (PExp)
protein for functional studies. Brand: TrueClones® Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Product Details	
Vector Backbone: pCMV6-XL5 Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studing the code of	Purpose:	Untagged full-length cDNA clone from Human SPDYE2L is ideal for over-expression of native protein for functional studies.
Promoter: Enhanced CMV Promoter, T7 Promoter Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Brand:	TrueClones®
Bacterial Resistance: Ampicillin Expression Type: Transient Characteristics: • These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies.	Vector Backbone:	pCMV6-XL5
Expression Type: Transient These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. These cDNA clones are ideal for over-expression of native proteins for functional studies.	Promoter:	Enhanced CMV Promoter, T7 Promoter
 Characteristics: These cDNA clones are isolated from full-length cDNA libraries and usually contain the coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. These cDNA clones are ideal for over-expression of native proteins for functional studies. 	Bacterial Resistance:	Ampicillin
 coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. These cDNA clones are ideal for over-expression of native proteins for functional studients. 	Expression Type:	Transient
 Every lot of primer is tested to provide clean sequencing of cDNA clones. 	Characteristics:	coding sequence as well as the untranslated regions (UTRs) of the mRNA transcript appropriate to the library from which they were isolated. • These cDNA clones are ideal for over-expression of native proteins for functional studies. Provided as 10 µg transfection-ready plasmids.
Purification: The DNAs were purified using PowerPrep HP Plasmid isolation kits for transfection ready plasmids.	Purification:	The DNAs were purified using PowerPrep HP Plasmid isolation kits for transfection ready plasmids.
Sequencing Primer: VP1.5 (forward) 5'GGACTTTCCAAAATGTCG 3', XL39 (reverse) 5'ATTAGGACAAGGCTGG	Sequencing Primer:	VP1.5 (forward) 5'GGACTTTCCAAAATGTCG 3', XL39 (reverse) 5'ATTAGGACAAGGCTGGTGGG

Product Details

Product Details	
	3'
Components:	 The cDNA clone is shipped in a 2-D bar-coded Matrix tube as dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Target Details	
Gene:	SPDYE2B
Alternative Name:	SPDYE2L (SPDYE2B Products)
NCBI Accession:	NM_001166339, NP_001159811
Application Details	
Restrictions:	For Research Use only
Handling	
Format:	Lyophilized
Storage:	RT,-20 °C
Storage Comment:	The lyophilized plasmid is stable for up to one year when stored at ambient temperature. Following dissolution in 100 μ L dH2O, store at -20 °C. Lyophilized primers are stable for up to one year when stored at ambient temperature. Following dissolution in 10 μ L dH2O, store at -20 °C.
Expiry Date:	12 months
Publications	
Product cited in:	Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (