-online.com **QENOMICS**



Datasheet for ABIN4677937

Mouse NME9 cDNA Clone in Bacterial Expression Vector (His-GST)

Quantity: 500 ng Gene: NME9 Species: Mouse Fusion tag: HIs-GST Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Mouse Txndc6 with HIs-GST Insert Length: 792 bp Vector Backbone: pPB-HIs-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 9-Xhol Fusion tag: Dual N-terminal tag, GX Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5-CACGTTTGGTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details Gene:	Overview	
Species: Mouse Fusion tag: His-GST Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Mouse Txndc6 with His-GST Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5-CACGTTTGGTGGTGGCGAC3, T7 terminator primer: 5-GCTAGTTATTGCTCAGCGG-3*	Quantity:	500 ng
Fusion tag: His-GST Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Mouse Txndc6 with His-GST Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag. Dual N-terminal tag. 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Gene:	NME9
Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Mouse Txndc6 with His-GST Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	Species:	Mouse
Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Mouse Txndc6 with His-GST Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Fusion tag:	His-GST
Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Mouse Txndc6 with His-GST Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Insert:	cDNA
Product Details Purpose: Bacterial expression of Mouse Txndc6 with His-GST Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Vector:	Bacterial Expression Vector
Purpose: Bacterial expression of Mouse Txndc6 with His-GST Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Application:	Cloning (Clon)
Insert Length: 792 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Product Details	
Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Purpose:	Bacterial expression of Mouse Txndc6 with His-GST
Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Insert Length:	792 bp
Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Vector Backbone:	pPB-His-GST
Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Promoter:	T7 Promoter
Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Bacterial Resistance:	Kanamycin
Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Expression Type:	Transient
which is cleavable with TEV (Size 27.9 kDa) Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details	Specificity:	5-Nhel and 3-Xhol
Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3' Target Details		Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein
GCTAGTTATTGCTCAGCGG-3' Target Details		which is cleavable with TEV (Size 27.9 kDa)
Target Details	Sequencing Primer:	GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-
		GCTAGTTATTGCTCAGCGG-3'
Gene: NME9	Target Details	
	Gene:	NME9

Target Details Txndc6 (NME9 Products) Alternative Name: NCBI Accession: NM_001165957 **Application Details Application Notes:** The pPB vectors are low-medium copy number vectors in which the gene expression is driven by the strong T7 promoter. Below are some basic guidelines for using the pPB vectors for protein production: 1. The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the host E. coli cell has a source of T7 RNA polymerase. 2. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl β-D-1thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM. 3. The ideal concentration of IPTG must be determined empirically for each recombinant protein/cell-line. Similarly, the length of time and temperature for induction provide other variables that need to be optimized on a case-to-case basis. 4. For toxic proteins, it is recommended to go for shorter induction time and also to try and suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins. 5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the purification. Restrictions: For Research Use only Handling

Format:	Liquid
Buffer:	10 mM Tris-HCI, 1 mM EDTA, pH 8.0
Storage:	-20 °C
Storage Comment:	1 year when stored at -20° C or lower in a non-frost free freezer.
Expiry Date:	12 months
5 1 11 11	

Publications

Product cited in:

Johnson, Drugan, Miller, Evans: "38" in:, Vol. 1363, Issue Nucleic acids research, pp. 28-39, (

1991)