-online.com **genomics**



Datasheet for ABIN4698196

Human BRAF cDNA Clone in Bacterial Expression Vector (His-MBP)

Quantity: 500 ng Gene: BRAF Species: Human Fusion tag: His-MBP Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Species: Human Fusion tag: His-MBP Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Fusion tag: His-MBP Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Purpose: Bacterial expression of Human BRAF with His-MBP Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Insert Length: 2299 bp Vector Backbone: pPB-His-MBP Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Vector Backbone:pPB-His-MBPPromoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-XholFusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose B	
	inding Protein which is
cleavable with Thrombin (Size 43 kDa)	
Sequencing Primer: MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 termin	ator primer: 5'-
GCTAGTTATTGCTCAGCGG-3'	
Target Details	
Gene: BRAF	

Target Details BRAF (BRAF Products) Alternative Name: NCBI Accession: NM_004333 **Application Details Application Notes:** The pPB vectors are low-medium copy number vectors in which the gene expression is driven by the strong T7 promoter. Below are some basic guidelines for using the pPB vectors for protein production: 1. The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the host E. coli cell has a source of T7 RNA polymerase. 2. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl β-D-1thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM. 3. The ideal concentration of IPTG must be determined empirically for each recombinant protein/cell-line. Similarly, the length of time and temperature for induction provide other variables that need to be optimized on a case-to-case basis. 4. For toxic proteins, it is recommended to go for shorter induction time and also to try and suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins. 5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the purification. Restrictions: For Research Use only Handling

Format:	Liquid
Buffer:	10 mM Tris-HCI, 1 mM EDTA, pH 8.0
Storage:	-20 °C
Storage Comment:	1 year when stored at -20° C or lower in a non-frost free freezer.
Expiry Date:	12 months

Publications

Product cited in:

Johnson, Drugan, Miller, Evans: "38" in:, Vol. 1363, Issue Nucleic acids research, pp. 28-39, (

1991)