-online.com genomics

Datasheet for ABIN4706288 Human MORF4L1 cDNA Clone in Bacterial Expression Vector (His-MBP)

Overview

| Quantity: | 500 ng |
|--------------|-----------------------------|
| Gene: | MORF4L1 |
| Species: | Human |
| Fusion tag: | His-MBP |
| Insert: | cDNA |
| Vector: | Bacterial Expression Vector |
| Application: | Cloning (Clon) |

Product Details

| Purpose:Bacterial expression of Human MORF4L1 with His-MBPInsert Length:972 bpVector Backbone:pPB-His-MBPPromoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:S-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa)Sequencing Primer:MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T2 terminator primer: 5'- CGTAGTTATTGCTCAGCGG-3' | | |
|---|-----------------------|--|
| Vector Backbone:pPB-His-MBPPromoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa)Sequencing Primer:MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 terminator primer: 5'- | Purpose: | Bacterial expression of Human MORF4L1 with His-MBP |
| Promoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa)Sequencing Primer:MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 terminator primer: 5'- | Insert Length: | 972 bp |
| Bacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa)Sequencing Primer:MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 terminator primer: 5'- | Vector Backbone: | pPB-His-MBP |
| Expression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa)Sequencing Primer:MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 terminator primer: 5'- | Promoter: | T7 Promoter |
| Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa) Sequencing Primer: MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 terminator primer: 5'- | Bacterial Resistance: | Kanamycin |
| Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa)Sequencing Primer:MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 terminator primer: 5'- | Expression Type: | Transient |
| | Specificity: | Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is |
| | Sequencing Primer: | |

Target Details

Gene:

MORF4L1

Alternative Name:

MORF4L1 (MORF4L1 Products)

Application Details

| Application Notes: | The pPB vectors are low-medium copy number vectors in which the gene expression is driven |
|--------------------|--|
| | by the strong T7 promoter. |
| | Below are some basic guidelines for using the pPB vectors for protein production: |
| | 1. The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the |
| | host E. coli cell has a source of T7 RNA polymerase. |
| | 2. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl β -D-1- |
| | thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM. |
| | 3. The ideal concentration of IPTG must be determined empirically for each recombinant |
| | protein/cell-line. Similarly, the length of time and temperature for induction provide other |
| | variables that need to be optimized on a case-to-case basis. |
| | 4. For toxic proteins, it is recommended to go for shorter induction time and also to try and |
| | suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS |
| | plasmid. Please note that special cell-lines are also available in the market that cater to |
| | expression of toxic proteins. |
| | 5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze |
| | the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the |
| | purification. |
| Restrictions: | For Research Use only |

Handling

| Format: | Liquid |
|-------------------|--|
| Buffer: | 10 mM Tris-HCI, 1 mM EDTA, pH 8.0 |
| Storage: | -20 °C |
| Storage Comment: | 1 year when stored at -20° C or lower in a non-frost free freezer. |
| Expiry Date: | 12 months |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (|
| | 1991) |