# -online.com genomics

## Datasheet for ABIN4709023 Human PRKCB cDNA Clone in Bacterial Expression Vector (His-MBP)

#### Overview

Quantity:	500 ng
Gene:	PKC beta (PRKCB)
Species:	Human
Fusion tag:	His-MBP
Insert:	cDNA
Vector:	Bacterial Expression Vector
Application:	Cloning (Clon)

#### Product Details

Vector Backbone:pPB-His-MBPPromoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein wh	)se:	Bacterial expression of Human PRKCB with His-MBP
Promoter:T7 PromoterBacterial Resistance:KanamycinExpression Type:TransientSpecificity:5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein wh	: Length:	2022 bp
Bacterial Resistance: Kanamycin   Expression Type: Transient   Specificity: 5-Nhel and 3-Xhol   Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein wh	or Backbone:	pPB-His-MBP
Expression Type: Transient   Specificity: 5-Nhel and 3-Xhol   Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein wh	oter:	T7 Promoter
Specificity: 5-Nhel and 3-Xhol   Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein wh	erial Resistance:	Kanamycin
Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein wh	ession Type:	Transient
cleavable with Thrombin (Size 43 kDa)	ficity:	5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Maltose Binding Protein which is cleavable with Thrombin (Size 43 kDa)
Sequencing Primer: MBP Forward primer: 5'-CGCAGATGTCCGCTTTCTGG-3', T7 terminator primer: 5'- GCTAGTTATTGCTCAGCGG-3'	encing Primer:	

### Target Details

Gene:

PKC beta (PRKCB)

Order at www.genomics-online.com

USA & Canada: +1 877 302 8632 | support@antibodies-online.com

Page 1/2 | Product datasheet for ABIN4709023 | 10/08/2023 | Copyright antibodies-online. All rights reserved.

Alternative Name:	PRKCB (PRKCB Products)
Application Details	
Application Notes:	The pPB vectors are low-medium copy number vectors in which the gene expression is driven
	by the strong T7 promoter.
	Below are some basic guidelines for using the pPB vectors for protein production:
	1. The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the
	host E. coli cell has a source of T7 RNA polymerase.
	2. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl $\beta$ -D-1-
	thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM.
	3. The ideal concentration of IPTG must be determined empirically for each recombinant
	protein/cell-line. Similarly, the length of time and temperature for induction provide other
	variables that need to be optimized on a case-to-case basis.
	4. For toxic proteins, it is recommended to go for shorter induction time and also to try and
	suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS
	plasmid. Please note that special cell-lines are also available in the market that cater to
	expression of toxic proteins.
	5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze
	the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the
	purification.
Restrictions:	For Research Use only
Handling	

Target Details

Format:	Liquid
Buffer:	10 mM Tris-HCl, 1 mM EDTA, pH 8.0
Storage:	-20 °C
Storage Comment:	1 year when stored at -20° C or lower in a non-frost free freezer.
Expiry Date:	12 months
Publications	
Product cited in:	Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (
	1991)