-online.com genomics

Datasheet for ABIN4827127 Human BAMBI cDNA Clone in Bacterial Expression Vector (His-GST)

Overview

Quantity:	500 ng
Gene:	BAMBI
Species:	Human
Fusion tag:	His-GST
Insert:	cDNA
Vector:	Bacterial Expression Vector
Application:	Cloning (Clon)

Product Details

Purpose:	Bacterial expression of Human BAMBI with His-GST
Insert Length:	783 bp
Vector Backbone:	pPB-His-GST
Promoter:	T7 Promoter
Bacterial Resistance:	Kanamycin
Expression Type:	Transient
Specificity:	5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Protein which is cleavable with TEV (Size 27.9 kDa)
Sequencing Primer:	GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'- GCTAGTTATTGCTCAGCGG-3'
T 10 1 1	

Target Details

0	~	201
G	ei	IE.

Alternative Name:

BAMBI (BAMBI Products)

Application Details

 4. For toxic proteins, it is recommended to go for shorter induction time and also to try and suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins. 5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the purification.
suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins. 5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze
suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins.
suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to
suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS
4. For toxic proteins, it is recommended to go for shorter induction time and also to try and
4 Endersia and in this account of the second state of the sheet of the structure in the second state of th
variables that need to be optimized on a case-to-case basis.
protein/cell-line. Similarly, the length of time and temperature for induction provide other
3. The ideal concentration of IPTG must be determined empirically for each recombinant
thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM.
2. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl β -D-1-
host E. coli cell has a source of T7 RNA polymerase.
1. The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the
Below are some basic guidelines for using the pPB vectors for protein production:
by the strong T7 promoter.
The pPB vectors are low-medium copy number vectors in which the gene expression is driven

Handling

Format:	Liquid
Buffer:	10 mM Tris-HCI, 1 mM EDTA, pH 8.0
Storage:	-20 °C
Storage Comment:	1 year when stored at -20° C or lower in a non-frost free freezer.
Expiry Date:	12 months
Publications	
Product cited in:	Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (
	1991)