-online.com **QENOMICS**



Datasheet for ABIN4831740

Human FPR3 cDNA Clone in Bacterial Expression Vector (His-GST)

Quantity: 500 ng Gene: FPR3 Species: Human Fusion tag: His-GST Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Species: Human Fusion tag: His-GST Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Fusion tag: His-GST Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Insert: cDNA Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Vector: Bacterial Expression Vector Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Application: Cloning (Clon) Product Details Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Purpose: Bacterial expression of Human FPR3 with His-GST Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Insert Length: 1062 bp Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Vector Backbone: pPB-His-GST Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Promoter: T7 Promoter Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Bacterial Resistance: Kanamycin Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Expression Type: Transient Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathion which is cleavable with TEV (Size 27.9 kDa)	
which is cleavable with TEV (Size 27.9 kDa)	
	e-S-Transferase Protein
Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminate	or primer: 5'-
GCTAGTTATTGCTCAGCGG-3'	
Target Details	
Gene: FPR3	

FPR3 (FPR3 Products)

Application Details

Application Notes:

The pPB vectors are low-medium copy number vectors in which the gene expression is driven by the strong T7 promoter.

Below are some basic guidelines for using the pPB vectors for protein production:

- 1. The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the host E. coli cell has a source of T7 RNA polymerase.
- 2. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl β-D-1thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM.
- 3. The ideal concentration of IPTG must be determined empirically for each recombinant protein/cell-line. Similarly, the length of time and temperature for induction provide other variables that need to be optimized on a case-to-case basis.
- 4. For toxic proteins, it is recommended to go for shorter induction time and also to try and suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins.
- 5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the purification.

Restrictions:

For Research Use only

Handling

Format:	Liquid
Buffer:	10 mM Tris-HCI, 1 mM EDTA, pH 8.0
Storage:	-20 °C
Storage Comment:	1 year when stored at -20° C or lower in a non-frost free freezer.
Expiry Date:	12 months

Publications

Product cited in:

Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)