## -online.com **QENOMICS**



Datasheet for ABIN4842667

## **Human TXN cDNA Clone in Bacterial Expression Vector (His-GST)**

Quantity:     500 ng       Gene:     TXN       Species:     Human       Fusion tag:     His-GST       Insert:     cDNA       Vector:     Bacterial Expression Vector       Application:     Cloning (Clon)       Product Details       Purpose:     Bacterial expression of Human TXN with His-GST       Insert Length:     318 bp       Vector Backbone:     pPB-His-GST       Promoter:     T7 Promoter       Bacterial Resistance:     Kanamycin       Expression Type:     Transient       Specificity:     5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoters and the properties of the propert	
Species: Human  Fusion tag: His-GST  Insert: cDNA  Vector: Bacterial Expression Vector  Application: Cloning (Clon)  Product Details  Purpose: Bacterial expression of Human TXN with His-GST  Insert Length: 318 bp  Vector Backbone: pPB-His-GST  Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Powhich is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Fusion tag: His-GST  Insert: cDNA  Vector: Bacterial Expression Vector  Application: Cloning (Clon)  Product Details  Purpose: Bacterial expression of Human TXN with His-GST  Insert Length: 318 bp  Vector Backbone: pPB-His-GST  Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoters: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Insert: cDNA  Vector: Bacterial Expression Vector  Application: Cloning (Clon)  Product Details  Purpose: Bacterial expression of Human TXN with His-GST  Insert Length: 318 bp  Vector Backbone: pPB-His-GST  Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoters: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Vector:  Application:  Cloning (Clon)  Product Details  Purpose:  Bacterial expression of Human TXN with His-GST  Insert Length:  318 bp  Vector Backbone:  pPB-His-GST  Promoter:  T7 Promoter  Bacterial Resistance:  Kanamycin  Expression Type:  Transient  Specificity:  5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoter is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer:  GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Application: Cloning (Clon)  Product Details  Purpose: Bacterial expression of Human TXN with His-GST  Insert Length: 318 bp  Vector Backbone: pPB-His-GST  Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoter is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Product Details  Purpose: Bacterial expression of Human TXN with His-GST  Insert Length: 318 bp  Vector Backbone: pPB-His-GST  Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Production of the Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Purpose: Bacterial expression of Human TXN with His-GST  Insert Length: 318 bp  Vector Backbone: pPB-His-GST  Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoter which is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCAGTTATTGCTCAGCGG-3'	
Insert Length: 318 bp  Vector Backbone: pPB-His-GST  Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoter which is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Vector Backbone:       pPB-His-GST         Promoter:       T7 Promoter         Bacterial Resistance:       Kanamycin         Expression Type:       Transient         Specificity:       5-Nhel and 3-Xhol         Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promotion in the property of the property	
Promoter: T7 Promoter  Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Promoter which is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Bacterial Resistance: Kanamycin  Expression Type: Transient  Specificity: 5-Nhel and 3-Xhol Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Province in the sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Expression Type:  Transient  Specificity:  5-Nhel and 3-Xhol  Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Province which is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer:  GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Specificity:  5-Nhel and 3-Xhol  Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Province of the Which is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer:  GST Forward primer: 5'-CACGTTTGGTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Fusion tag: Dual N-terminal tag, 6X Histidine followed by Glutathione-S-Transferase Provided Primer:  Which is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer:  GST Forward primer: 5'-CACGTTTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
which is cleavable with TEV (Size 27.9 kDa)  Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	
Sequencing Primer: GST Forward primer: 5'-CACGTTTGGTGGCGAC3', T7 terminator primer: 5'-GCTAGTTATTGCTCAGCGG-3'	rotein
GCTAGTTATTGCTCAGCGG-3'	
Target Details	
Gene: TXN	

## **Target Details** TXN (TXN Products) Alternative Name: NCBI Accession: NM\_003329 **Application Details Application Notes:** The pPB vectors are low-medium copy number vectors in which the gene expression is driven by the strong T7 promoter. Below are some basic guidelines for using the pPB vectors for protein production: 1. The pPB vectors are designed to be used with E. coli strains that are DE3 lysogens i.e. the host E. coli cell has a source of T7 RNA polymerase. 2. Recombinant protein induction is usually done at OD600 of 0.6-1.2 using Isopropyl β-D-1thiogalactopyranoside (IPTG) at a final concentration of 0.05 -1mM. 3. The ideal concentration of IPTG must be determined empirically for each recombinant protein/cell-line. Similarly, the length of time and temperature for induction provide other variables that need to be optimized on a case-to-case basis. 4. For toxic proteins, it is recommended to go for shorter induction time and also to try and suppress basal recombinant gene expression through (a) addition of glucose or use of pLysS plasmid. Please note that special cell-lines are also available in the market that cater to expression of toxic proteins. 5. Once grown for the desired length of time, harvest cells by centrifugation and either freeze the cells at -80°C (as such or after re-suspending in the desired buffer) or proceed with the purification. Restrictions: For Research Use only Handling

Format:	Liquid
Buffer:	10 mM Tris-HCl, 1 mM EDTA, pH 8.0
Storage:	-20 °C
Storage Comment:	1 year when stored at -20° C or lower in a non-frost free freezer.
Expiry Date:	12 months

## **Publications**

Product cited in:

Johnson, Drugan, Miller, Evans: "38" in:, Vol. 1363, Issue Nucleic acids research, pp. 28-39, (

1991)