-online.com genomics





Human FAM183A ORF Clone in Mammalian Expression Vector (DYKDDDDK Tag)

Overview	
Quantity:	10 μg
Gene:	FAM183A
Species:	Human
Fusion tag:	DYKDDDDK Tag
Insert:	ORF
Vector:	Mammalian Expression Vector
Application:	Protein Expression (PExp)
Product Details	
Purpose:	Expression/transfection ready cDNA ORF clone of Human FAM183A with C terminal
	DYKDDDDK tag is ideal for express proteins in E.coli & mammalian cells.
Brand:	GenEZ™
Insert Length:	405 bp
Vector Backbone:	pcDNA3.1+C-(K)-DYK
Promoter:	CMV Promoter
Selectable Marker:	Neomycin
Bacterial Resistance:	Ampicillin
Expression Type:	Transient, Stable
Sequence:	ATGGCGGGAC ACCCAAAAGA GAAGGTGATT CCAGATGAGG TCCATCAGAA CCAGATCTTG
	CGGGAACTGT ACCTCAAGGA GCTACGAACC CAGAAACTCT ACACGCAGTA TCACGTGAAT
	CCCCTCCGCA AGATTCATAC AGTCACCAGG AAGCCCATGT CTTGGCATGA TAACCTGGAG
	GAACCTGCAG ATGCCAGGTT TCTGAATCTC ATTCACCATG CTGCCCAGGG ACCAAGGAAG

Product Details

	AAGTACCCAG AGACACAGAC TGAAAACCAG GAAGTTGGAT GGGACTTAGA GCCCTTGATC
	AACCCAGAAC GCCATGACCG CAGGCTGAAT CACTTCAGGG TCTGCAGTGA CATCACTCTG
	TACAAGGCTA AAACGTGGGG CTTAGGAGAT GATCACCACA AGTAG
Specificity:	ORF Insert Method: CloneEZ® Seamless cloning technology, recombination-based cloning technology
Characteristics:	Gene cDNA ORF clone sequences were retrieved from the NCBI Reference Sequence Database (RefSeq). These sequences represent the protein coding region of the gene cDNA ORF which is encoded by the open reading frame (ORF) sequence.
Sequencing Primer:	 Forward primer: 5'-TAATACGACTCACTATAGGG-3' Reverse primer: 5'-CCTCGACTGTGCCTTCTA-3'
Grade:	End-sequenced
Components:	The GenEZ ORF clone is delivered as 10 µg of lyophilized plasmid DNA in a vial.
Target Details	
Gene:	FAM183A
Alternative Name:	FAM183A (FAM183A Products)
Gene ID:	440585
NCBI Accession:	NM_001101376
Application Details	
Restrictions:	For Research Use only
Handling	
Format:	Lyophilized
Storage:	RT/-20 °C
Storage Comment:	 Keep the vial sealed and store at -20°C for long-term storage. Before use, centrifuge the vial at 6,000 g x g for 1 minute at 4°C. Open the lid and add 100 µl (or other volume depending on your desired final concentration) of distilled water (or TE buffer) to dissolve the DNA. If necessary, heat the solution at 50°C for 15 minutes to dissolve the DNA. Close the lid and vortex the vial for 1 minute. Aliquot the dissolved plasmid DNA and store in small aliquots at -20°C.

Handling	
Expiry Date:	12 months
Publications	
Product cited in:	Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (
	1991)