-online.com **Genomics**





Human TAS2R16 CRISPR gRNA in Lenti Particles

| Overview | |
|-----------------------|-----------------------------------------------------------------------------------------------------------------|
| Quantity: | 3 x 300 μL |
| Gene: | TAS2R16 |
| Species: | Human |
| Insert: | gRNA |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Set of 3 gRNA against TAS2R16 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (Cas9 required separately) |
| Vector Backbone: | pLenti-U6-sgRNA-PGK-Neo |
| Promoter: | U6 Promoter, PGK Promoter |
| Selectable Marker: | Neomycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |
| Sequence: | Sequence available upon placing order |
| Specificity: | GRNAs are designed for use with Cas9 Nuclease only. |
| Components: | Lentiviral particles with a set of 3 gRNAs (3 x 300 μ L) covering different sequences of TAS2R16 |
| Target Details | |
| Gene: | TAS2R16 |
| Alternative Name: | TAS2R16 (TAS2R16 Products) |
| NCBI Accession: | NM_016945 |

Application Details

| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing |
|--------------------|--------------------------------------------------------------------------------------------------|
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Viral Particles |
| Storage: | -80 °C |
| Expiry Date: | 12 months |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (|