-online.com Jenomics

Datasheet for ABIN5166562 Human TMEM129 CRISPR gRNA in Lenti Particles

| Overview | |
|---------------------------------|---|
| Quantity: | 300 µL |
| Gene: | TMEM129 |
| Species: | Human |
| Insert: | gRNA |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Individual gRNA against TMEM129 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (Cas9 required separately) |
| Vector Backbone: | |
| | pLenti-U6-sgRNA-PGK-Neo |
| Promoter: | pLenti-U6-sgRNA-PGK-Neo U6 Promoter, PGK Promoter |
| | |
| Promoter: | U6 Promoter, PGK Promoter |
| Promoter: Selectable Marker: | U6 Promoter, PGK Promoter Neomycin |

| Specificity: | GRNAs are designed for use with Cas9 Nuclease only. |
|--------------|---|
| Components: | Lentiviral particles with an individual gRNA (300 μ L) for a specific sequence of TMEM129 |

Target Details

| Gene: | TMEM129 |
|-------------------|----------------------------|
| Alternative Name: | TMEM129 (TMEM129 Products) |
| NCBI Accession: | NM_001127266 |

Order at www.genomics-online.com USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN5166562 | 09/12/2023 | Copyright antibodies-online. All rights reserved.

| Application Details | |
|---------------------|---|
| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Viral Particles |
| Storage: | -80 °C |
| Expiry Date: | 12 months |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991) |