-online.com Genomics



| Overview | |
|-----------------------|--|
| Quantity: | 3 x 300 µL |
| Gene: | USF3 (KIAA2018) |
| Species: | Human |
| Insert: | gRNA |
| Vector: | Lentiviral Vector |
| Application: | Genome Editing with Engineered Nucleases (GEEN), Protein Expression (PExp) |
| Product Details | |
| Purpose: | Set of 3 gRNA against KIAA2018 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (Cas9 required separately) |
| Vector Backbone: | pLenti-U6-sgRNA-PGK-Neo |
| Promoter: | U6 Promoter, PGK Promoter |
| Selectable Marker: | Neomycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |
| Sequence: | Sequence available upon placing order |
| Specificity: | GRNAs are designed for use with Cas9 Nuclease only. |
| Components: | Lentiviral particles with a set of 3 gRNAs (3 x 300 $\mu L)$ covering different sequences of KIAA2018 |

Target Details

| Gene: | USF3 (KIAA2018) |
|-------------------|-----------------|
| Alternative Name: | KIAA2018 |
| NCBI Accession: | NM_001009899 |

Order at www.genomics-online.com USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN5176409 | 09/12/2023 | Copyright antibodies-online. All rights reserved.

| Application Details | |
|---------------------|---|
| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Viral Particles |
| Storage: | -80 °C |
| Expiry Date: | 12 months |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991) |