## -online.com **Genomics**





## **Human BVES-AS1 CRISPR gRNA in Lenti Particles**

| Overview              |  |
|-----------------------|--|
| Quantity:             | 3 x 300 μL   |
| Gene:                 | C6orf112 (BVES-AS1)  |
| Species:              | Human  |
| Insert:               | gRNA   |
| Vector:               | Lentiviral Vector  |
| Application:          | Genome Editing with Engineered Nucleases (GEEN), Protein Expression (PExp)                                       |
| Product Details       |  |
| Purpose:              | Set of 3 gRNA against C6orf112 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (Cas9 required separately) |
| Vector Backbone:      | pLenti-U6-sgRNA-PGK-Neo  |
| Promoter:             | U6 Promoter, PGK Promoter  |
| Selectable Marker:    | Neomycin   |
| Bacterial Resistance: | Ampicillin   |
| Expression Type:      | Stable, Transient  |
| Sequence:             | Sequence available upon placing order  |
| Specificity:          | GRNAs are designed for use with Cas9 Nuclease only.  |
| Components:           | Lentiviral particles with a set of 3 gRNAs (3 x 300 $\mu$ L) covering different sequences of C6orf112            |
| Target Details        |  |
| Gene:                 | C6orf112 (BVES-AS1)  |
| Alternative Name:     | C6orf112   |

## **Application Details**

| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing                        |
|--------------------|--|
| Restrictions:      | For Research Use only  |
| Handling           |  |
| Format:            | Viral Particles  |
| Storage:           | -80 °C   |
| Expiry Date:       | 12 months  |
| Publications       |  |
| Product cited in:  | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, ( |