## -online.com **Genomics**





## Human CDK3 CRISPR gRNA + Cas9 in Lenti Particles

| Overview              |   |
|-----------------------|---|
| Quantity:             | 3 x 300 μL  |
| Gene:                 | CDK3  |
| Species:              | Human   |
| Insert:               | gRNA + Cas9   |
| Vector:               | Lentiviral Vector   |
| Application:          | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN)                        |
| Product Details       |   |
| Purpose:              | Set of 3 gRNA against CDK3 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and      |
|                       | Cas9 in a single vector)  |
| Vector Backbone:      | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro   |
| Promoter:             | U6 Promoter, SFFV Promoter  |
| Selectable Marker:    | Puromycin   |
| Bacterial Resistance: | Ampicillin  |
| Expression Type:      | Stable, Transient   |
| Sequence:             | Sequence available upon placing order   |
| Specificity:          | GRNAs are designed for use with Cas9 Nuclease only.   |
|                       | Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority     |
|                       | of cells, except ES cells or iPS cells.   |
| Sequencing Primer:    | U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA-3'  |
| Components:           | Lentiviral particles with a set of 3 gRNAs (3 x 300 $\mu$ L) covering different sequences of CDK3 |

## **Target Details** CDK3 Gene: Alternative Name: CDK3 (CDK3 Products) NCBI Accession: NM\_001258 **Application Details** Application Notes: Recommended for quality control: Restriction Enzyme Digest and Sequencing Restrictions: For Research Use only Handling Viral Particles Format: -80 °C Storage: Expiry Date: 12 months

## **Publications**

Product cited in: Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (
1991)