## -online.com **Genomics**

Specificity:

Sequencing Primer:

Components:

## Datasheet for ABIN5221000 Human IFT172 CRISPR gRNA + Cas9 in Lenti Particles

Overview	
Quantity:	300 µL
Gene:	IFT172
Species:	Human
Insert:	gRNA + Cas9
Vector:	Lentiviral Vector
Application:	Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN)
Product Details	
Purpose:	Individual gRNA against IFT172 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector)
Vector Backbone:	pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro
Promoter:	U6 Promoter, SFFV Promoter
Selectable Marker:	Puromycin
Bacterial Resistance:	Ampicillin
Expression Type:	Stable, Transient
Sequence:	Sequence available upon placing order

U6 Forward Primer: 5'--TACGTCCAAGGTCGGGCAGGAAGA--3'

GRNAs are designed for use with Cas9 Nuclease only.

of cells, except ES cells or iPS cells.

Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority

Lentiviral particles with an individual gRNA (300  $\mu$ L) for a specific sequence of IFT172

Target Details	
Gene:	IFT172
Alternative Name:	IFT172 (IFT172 Products)
NCBI Accession:	NM_015662
Application Details	
Application Notes:	Recommended for quality control: Restriction Enzyme Digest and Sequencing
Restrictions:	For Research Use only
Handling	
Format:	Viral Particles
Storage:	-80 °C
Expiry Date:	12 months
Publications	
Product cited in:	Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, ( 1991)