

Datasheet for ABIN5240015

Mouse MAPK8 CRISPR gRNA + Cas9 in Lenti Particles

Overview

| | |
|--------------|--|
| Quantity: | 300 µL |
| Gene: | JNK (MAPK8) |
| Species: | Mouse |
| Insert: | gRNA + Cas9 |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PEXP), Genome Editing with Engineered Nucleases (GEEN) |

Product Details

| | |
|-----------------------|--|
| Purpose: | Individual gRNA against Mapk8 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector) |
| Vector Backbone: | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro |
| Promoter: | U6 Promoter, SFFV Promoter |
| Selectable Marker: | Puromycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |
| Sequence: | Sequence available upon placing order |
| Specificity: | GRNAs are designed for use with Cas9 Nuclease only. Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority of cells, except ES cells or iPS cells. |
| Sequencing Primer: | U6 Forward Primer: 5'--TACGTCCAAGGTCGGGCAGGAAGA--3' |
| Components: | Lentiviral particles with an individual gRNA (300 µL) for a specific sequence of Mapk8 |

Target Details

Gene: JNK (MAPK8)

Alternative Name: Mapk8 ([MAPK8 Products](#))

NCBI Accession: [NM_016700](#)

Application Details

Application Notes: Recommended for quality control: Restriction Enzyme Digest and Sequencing

Restrictions: For Research Use only

Handling

Format: Viral Particles

Storage: -80 °C

Expiry Date: 12 months

Publications

Product cited in: Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)