-online.com Jenomics

Datasheet for ABIN5251212 Human TP53BP2 CRISPR gRNA + Cas9 in Lenti Particles

3 x 300 µL
TP53BP2
Human
gRNA + Cas9
Lentiviral Vector
Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN)
Set of 3 gRNA against TP53BP2 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA
and Cas9 in a single vector)
pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro
U6 Promoter, SFFV Promoter
Puromycin
Ampicillin
Stable, Transient
Sequence available upon placing order
GRNAs are designed for use with Cas9 Nuclease only.
Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority
of cells, except ES cells or iPS cells.
U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA3'

Components:

Lentiviral particles with a set of 3 gRNAs (3 x 300 µL) covering different sequences of TP53BP2

Target Details	
Gene:	TP53BP2
Alternative Name:	TP53BP2 (TP53BP2 Products)
NCBI Accession:	NM_001031685
Application Details	
Application Notes:	Recommended for quality control: Restriction Enzyme Digest and Sequencing
Restrictions:	For Research Use only
Handling	
Format:	Viral Particles
Storage:	-80 °C
Expiry Date:	12 months
Publications	
Product cited in:	Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)