genomics -online.com





Human KDM5D CRISPR gRNA + Cas9 in Lenti Particles

Overview	
Quantity:	3 x 300 μL
Gene:	KDM5D
Species:	Human
Insert:	gRNA + Cas9
Vector:	Lentiviral Vector
Application:	Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN)
Product Details	
Purpose:	Set of 3 gRNA against KDM5D in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and
	Cas9 in a single vector)
Vector Backbone:	pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro
Promoter:	U6 Promoter, SFFV Promoter
Selectable Marker:	Puromycin
Bacterial Resistance:	Ampicillin
Expression Type:	Stable, Transient
Sequence:	Sequence available upon placing order
Specificity:	GRNAs are designed for use with Cas9 Nuclease only.
	Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority
	of cells, except ES cells or iPS cells.
Sequencing Primer:	U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA3'
Components:	Lentiviral particles with a set of 3 gRNAs (3 x 300 μ L) covering different sequences of KDM5D

Target Details

Gene:	KDM5D
Alternative Name:	KDM5D (KDM5D Products)
NCBI Accession:	NM_001146705

Application Notes:	Recommended for quality control: Restriction Enzyme Digest and Sequencing
Restrictions:	For Research Use only
Handling	
Format:	Viral Particles
Storage:	-80 °C
Expiry Date:	12 months
Publications	
Product cited in:	Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39

Product cited in: Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (
1991)