-online.com Genomics

Datasheet for ABIN5257264 Mouse GTF3C1 CRISPR gRNA + Cas9 in Lenti Particles

| Overview | |
|-----------------------|--|
| Quantity: | 3 x 300 μL |
| Gene: | GTF3C1 |
| Species: | Mouse |
| Insert: | gRNA + Cas9 |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Set of 3 gRNA against Gtf3c1 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and |
| | Cas9 in a single vector) |
| Vector Backbone: | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro |
| Promoter: | U6 Promoter, SFFV Promoter |
| Selectable Marker: | Puromycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |
| Sequence: | Sequence available upon placing order |

Specificity:GRNAs are designed for use with Cas9 Nuclease only.
Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority
of cells, except ES cells or iPS cells.Sequencing Primer:U6 Forward Primer: 5'--TACGTCCAAGGTCGGGCAGGAAGA--3'Components:Lentiviral particles with a set of 3 gRNAs (3 x 300 µL) covering different sequences of Gtf3c1

| Target Details | |
|---------------------|--|
| Gene: | GTF3C1 |
| Alternative Name: | Gtf3c1 (GTF3C1 Products) |
| NCBI Accession: | NM_207239 |
| Application Details | |
| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Viral Particles |
| Storage: | -80 °C |
| Expiry Date: | 12 months |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (|
| | 1991) |