-online.com genomics

Datasheet for ABIN5257609 Human GLDC CRISPR gRNA + Cas9 in Lenti Particles

| Overview | |
|-----------------------|--|
| Quantity: | 300 µL |
| Gene: | GLDC |
| Species: | Human |
| Insert: | gRNA + Cas9 |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Individual gRNA against GLDC in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and |
| | Cas9 in a single vector) |
| Vector Backbone: | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro |
| Promoter: | U6 Promoter, SFFV Promoter |
| Selectable Marker: | Puromycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |

U6 Forward Primer: 5'--TACGTCCAAGGTCGGGCAGGAAGA--3'

Sequence available upon placing order

of cells, except ES cells or iPS cells.

GRNAs are designed for use with Cas9 Nuclease only.

Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority

Lentiviral particles with an individual gRNA (300 μ L) for a specific sequence of GLDC

Sequence:

Specificity:

Sequencing Primer:

Components:

| Target Details | |
|---------------------|--|
| Gene: | GLDC |
| Alternative Name: | GLDC (GLDC Products) |
| NCBI Accession: | NM_000170 |
| Application Details | |
| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Viral Particles |
| Storage: | -80 °C |
| Expiry Date: | 12 months |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (|
| | 1991) |