## -online.com genomics

## Datasheet for ABIN5262789 Human GLS2 CRISPR gRNA + Cas9 in Lenti Particles

| Over | rview |
|------|-------|
| Ovei | VIEW  |

| Quantity:    | 300 µL   |
|--------------|--|
| Gene:        | GLS2   |
| Species:     | Human  |
| Insert:      | gRNA + Cas9  |
| Vector:      | Lentiviral Vector  |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |

## Product Details

| Purpose:              | Individual gRNA against GLS2 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector)   |
|-----------------------|---|
| Vector Backbone:      | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro   |
| Promoter:             | U6 Promoter, SFFV Promoter  |
| Selectable Marker:    | Puromycin   |
| Bacterial Resistance: | Ampicillin  |
| Expression Type:      | Stable, Transient   |
| Sequence:             | Sequence available upon placing order   |
| Specificity:          | GRNAs are designed for use with Cas9 Nuclease only.<br>Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority<br>of cells, except ES cells or iPS cells. |
| Sequencing Primer:    | U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA3'   |
| Components:           | Lentiviral particles with an individual gRNA (300 $\mu L)$ for a specific sequence of GLS2  |

| Target Details      |  |
|---------------------|--|
| Gene:               | GLS2   |
| Alternative Name:   | GLS2 (GLS2 Products)   |
| NCBI Accession:     | NM_013267  |
| Application Details |  |
| Application Notes:  | Recommended for quality control: Restriction Enzyme Digest and Sequencing                        |
| Restrictions:       | For Research Use only  |
| Handling            |  |
| Format:             | Viral Particles  |
| Storage:            | -80 °C   |
| Expiry Date:        | 12 months  |
| Publications        |  |
| Product cited in:   | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, ( |
|                     | 1991)  |