-online.com genomics

Datasheet for ABIN5271665 Human NAV3 CRISPR gRNA + Cas9 in Lenti Particles

Overview	

Quantity:	300 µL
Gene:	NAV3
Species:	Human
Insert:	gRNA + Cas9
Vector:	Lentiviral Vector
Application:	Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN)

Product Details

Purpose:	Individual gRNA against NAV3 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector)
Vector Backbone:	pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro
Promoter:	U6 Promoter, SFFV Promoter
Selectable Marker:	Puromycin
Bacterial Resistance:	Ampicillin
Expression Type:	Stable, Transient
Sequence:	Sequence available upon placing order
Specificity:	GRNAs are designed for use with Cas9 Nuclease only. Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority of cells, except ES cells or iPS cells.
Sequencing Primer:	U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA3'
Components:	Lentiviral particles with an individual gRNA (300 μ L) for a specific sequence of NAV3

NAV3
NAV3 (NAV3 Products)
NM_014903
Recommended for quality control: Restriction Enzyme Digest and Sequencing
For Research Use only
Viral Particles
-80 °C
12 months
Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (