-online.com Genomics

Specificity:

Sequencing Primer:

Components:

Datasheet for ABIN5275759 Human TMEM40 CRISPR gRNA + Cas9 in Lenti Particles

| Overview | |
|-----------------------|---|
| Quantity: | 300 µL |
| Gene: | TMEM40 |
| Species: | Human |
| Insert: | gRNA + Cas9 |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Individual gRNA against TMEM40 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector) |
| Vector Backbone: | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro |
| Promoter: | U6 Promoter, SFFV Promoter |
| Selectable Marker: | Puromycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |
| Sequence: | Sequence available upon placing order |

U6 Forward Primer: 5'--TACGTCCAAGGTCGGGCAGGAAGA--3'

GRNAs are designed for use with Cas9 Nuclease only.

of cells, except ES cells or iPS cells.

Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority

Lentiviral particles with an individual gRNA (300 μ L) for a specific sequence of TMEM40

| Target Details | |
|---------------------|---|
| Gene: | TMEM40 |
| Alternative Name: | TMEM40 (TMEM40 Products) |
| NCBI Accession: | NM_018306 |
| Application Details | |
| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing |
| Restrictions: | For Research Use only |
| Handling | |
| Format: | Viral Particles |
| Storage: | -80 °C |
| Expiry Date: | 12 months |
| Publications | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991) |