-online.com

Datasheet for ABIN5281415 Mouse OLFR146 CRISPR gRNA + Cas9 in Lenti Particles

| Overview | |
|-----------------------|---|
| Quantity: | 300 µL |
| Gene: | Olfr146 (OLFR146) |
| Species: | Mouse |
| Insert: | gRNA + Cas9 |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Individual gRNA against Olfr146 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector) |
| Vector Backbone: | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro |
| Promoter: | U6 Promoter, SFFV Promoter |
| Selectable Marker: | Puromycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |
| Sequence: | Sequence available upon placing order |
| Specificity: | GRNAs are designed for use with Cas9 Nuclease only. Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority of cells, except ES cells or iPS cells. |
| Sequencing Primer: | U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA3' |
| Components: | Lentiviral particles with an individual gRNA (300 μ L) for a specific sequence of Olfr146 |

| Target Details | | |
|---------------------|--|--|
| Gene: | Olfr146 (OLFR146) | |
| Alternative Name: | Olfr146 | |
| NCBI Accession: | NM_146747 | |
| Application Details | | |
| Application Notes: | Recommended for quality control: Restriction Enzyme Digest and Sequencing | |
| Restrictions: | For Research Use only | |
| Handling | | |
| Format: | Viral Particles | |
| Storage: | -80 °C | |
| Expiry Date: | 12 months | |
| Publications | | |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (| |
| | 1991) | |