genomics -online.com



Datasheet for ABIN5295382

Human TMEM133 CRISPR gRNA + Cas9 in Lenti Particles

| Overview | |
|-----------------------|--|
| Quantity: | 3 x 300 μL |
| Gene: | TMEM133 |
| Species: | Human |
| Insert: | gRNA + Cas9 |
| Vector: | Lentiviral Vector |
| Application: | Protein Expression (PExp), Genome Editing with Engineered Nucleases (GEEN) |
| Product Details | |
| Purpose: | Set of 3 gRNA against TMEM133 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector) |
| Vector Backbone: | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro |
| Promoter: | U6 Promoter, SFFV Promoter |
| Selectable Marker: | Puromycin |
| Bacterial Resistance: | Ampicillin |
| Expression Type: | Stable, Transient |
| Sequence: | Sequence available upon placing order |
| Specificity: | GRNAs are designed for use with Cas9 Nuclease only. Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority of cells, except ES cells or iPS cells. |
| Sequencing Primer: | U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA3' |
| Components: | Lentiviral particles with a set of 3 gRNAs (3 x 300 μ L) covering different sequences of TMEM133 |

Target Details TMEM133 Gene: Alternative Name: TMEM133 NCBI Accession: NM_032021 **Application Details** Application Notes: Recommended for quality control: Restriction Enzyme Digest and Sequencing Restrictions: For Research Use only Handling Viral Particles Format: -80 °C Storage:

Expiry Date:

Product cited in: Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (
1991)

12 months