## -online.com genomics





## Human PTPRZ2 CRISPR gRNA + Cas9 in Lenti Particles

| Overview              |  |
|-----------------------|--|
| Quantity:             | 300 μL   |
| Gene:                 | PTPRZ (PTPRZ2)   |
| Species:              | Human  |
| Insert:               | gRNA + Cas9  |
| Vector:               | Lentiviral Vector  |
| Application:          | Genome Editing with Engineered Nucleases (GEEN), Protein Expression (PExp)   |
| Product Details       |  |
| Purpose:              | Individual gRNA against PTPRZ2 in Lentiviral Particles with a Titer of >1x10e7 IU/mL. (sgRNA and Cas9 in a single vector)  |
| Vector Backbone:      | pLenti-U6-sgRNA-SFFV-Cas9-2A-Puro  |
| Promoter:             | U6 Promoter, SFFV Promoter   |
| Selectable Marker:    | Puromycin  |
| Bacterial Resistance: | Ampicillin   |
| Expression Type:      | Stable, Transient  |
| Sequence:             | Sequence available upon placing order  |
| Specificity:          | GRNAs are designed for use with Cas9 Nuclease only.  Cas9 Nuclease is under the control of the SFFV promoter which should work for a vast majority of cells, except ES cells or iPS cells. |
| Sequencing Primer:    | U6 Forward Primer: 5'TACGTCCAAGGTCGGGCAGGAAGA3'  |
| Components:           | Lentiviral particles with an individual gRNA (300 $\mu$ L) for a specific sequence of PTPRZ2   |

## Gene: PTPRZ (PTPRZ2) Alternative Name: PTPRZ2 (PTPRZ2 Products) Application Details Application Notes: Recommended for quality control: Restriction Enzyme Digest and Sequencing Restrictions: For Research Use only Handling

| Restrictions:     | For Research Use only  |
|-------------------|--|
| Handling          |  |
| Format:           | Viral Particles  |
| Storage:          | -80 °C   |
| Expiry Date:      | 12 months  |
| Publications      |  |
| Product cited in: | Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, ( |

1991)