## -online.com genomics



Overview	
Quantity:	10 µg
Insert:	Empty
Vector:	Lentiviral Vector
Application:	Cloning (Clon), Protein Extraction (PEx)
Product Details	
Purpose:	Clone your gene of interest into this Lentiviral Expression Vector, then co-transfect along with lentiviral packaging vectors.
Vector Backbone:	pSMPUW-Hygro
Promoter:	EFS Promoter, PGK Promoter
Selectable Marker:	Hygromycin
Bacterial Resistance:	Kanamycin
Expression Type:	Stable, Transient
Application Details	
Application Notes:	Optimal working dilution should be determined by the investigator.
Restrictions:	For Research Use only
Handling	
Format:	Lyophilized
Storage:	-20 °C

Publications	
--------------	--

Product cited in:

Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, ( 1991)