

Datasheet for ABIN5691569 pSMPUW-Puro Lentiviral Expression Vector

Overview

Quantity:	10 µg
Insert:	Empty
Vector:	Lentiviral Vector
Application:	Cloning (Clon), Protein Extraction (PEx)

Product Details

Purpose:	Clone your gene of interest into this Lentiviral Expression Vector, then co-transfect along with lentiviral packaging vectors.
Vector Backbone:	pSMPUW-Puro
Promoter:	PGK Promoter, EFS Promoter
Selectable Marker:	Puromycin
Bacterial Resistance:	Kanamycin
Expression Type:	Stable, Transient

Application Details

Application Notes:	Optimal working dilution should be determined by the investigator.
Restrictions:	For Research Use only

Handling

Format:	Lyophilized
Storage:	-20 °C

Publications

Product cited in: Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)