-online.com genomics

Datasheet for ABIN6391940 Human NPC2 ORF Clone in Lenti Particles (GFP tag)

Overview

Quantity:	200 μL
Gene:	NPC2
Species:	Human
Fusion tag:	GFP tag
Insert:	ORF
Vector:	Lentiviral Vector
Application:	Protein Expression (PExp)

Product Details

Vector Backbone: pLenti-C-mGFP Promoter: CMV Promoter Bacterial Resistance: Chloramphenicol Expression Type: Transient Specificity: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. Characteristics: • Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect & hard-to-transfect cells. • High transduction efficiency. • Convenience: Minimal need for optimization. • Safety: 3rd generation system with improved biosafety. • Pre-titered, ready-to-use • Titer guaranteed, 10^7 TU/mL • Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity		
Promoter: CMV Promoter Bacterial Resistance: Chloramphenicol Expression Type: Transient Specificity: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. Characteristics: Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect & hard-to-transfect cells. Characteristics: High transduction efficiency. Convenience: Minimal need for optimization. Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity Components: Lentiviral particles with guaranteed titer of >10^7 TU/mL	Purpose:	Lenti ORF particles, NPC2 (mGFP-tagged) - Human Niemann-Pick disease, type C2 (NPC2)
Bacterial Resistance: Chloramphenicol Expression Type: Transient Specificity: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. Characteristics: Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect & hard-to-transfect cells. High transduction efficiency. Convenience: Minimal need for optimization. Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity Components: Lentiviral particles with guaranteed titer of >10^7 TU/mL	Vector Backbone:	pLenti-C-mGFP
Expression Type:TransientSpecificity:This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.Characteristics:Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect & hard-to-transfect cells.High transduction efficiency.Convenience: Minimal need for optimization.Safety: 3rd generation system with improved biosafety.Pre-titered, ready-to-useTiter guaranteed, 10^7 TU/mLProvided in the proprietary Lenti Stabilizer Solution with 1 year infectivity	Promoter:	CMV Promoter
Specificity: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. Characteristics: • Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect & hard-to-transfect cells. • High transduction efficiency. • Convenience: Minimal need for optimization. • Safety: 3rd generation system with improved biosafety. • Pre-titered, ready-to-use • Titer guaranteed, 10^7 TU/mL • Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity Components: Lentiviral particles with guaranteed titer of >10^7 TU/mL	Bacterial Resistance:	Chloramphenicol
varies depending on the nature of the gene. Characteristics: Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect & hard-to-transfect cells. High transduction efficiency. Convenience: Minimal need for optimization. Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity Components:	Expression Type:	Transient
 Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect & hard-to-transfect cells. High transduction efficiency. Convenience: Minimal need for optimization. Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity 	Specificity:	This clone was engineered to express the complete ORF with an expression tag. Expression
hard-to-transfect cells. High transduction efficiency. Convenience: Minimal need for optimization. Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity Components: Lentiviral particles with guaranteed titer of >10^7 TU/mL		varies depending on the nature of the gene.
 High transduction efficiency. Convenience: Minimal need for optimization. Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity 	Characteristics:	• Broad cell spectrum: Lentivirus infect many cells, dividing & non-dividing, easy-to-transfect &
 Convenience: Minimal need for optimization. Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity 		hard-to-transfect cells.
 Safety: 3rd generation system with improved biosafety. Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity 		High transduction efficiency.
 Pre-titered, ready-to-use Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity 		Convenience: Minimal need for optimization.
 Titer guaranteed, 10^7 TU/mL Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity Components: Lentiviral particles with guaranteed titer of >10^7 TU/mL 		 Safety: 3rd generation system with improved biosafety.
 Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity Components: Lentiviral particles with guaranteed titer of >10^7 TU/mL 		Pre-titered, ready-to-use
Components: Lentiviral particles with guaranteed titer of >10^7 TU/mL		Titer guaranteed, 10^7 TU/mL
		Provided in the proprietary Lenti Stabilizer Solution with 1 year infectivity
Order at www.genomics-online.com	Components:	Lentiviral particles with guaranteed titer of >10^7 TU/mL

USA & Canada: +1 877 302 8632 | support@antibodies-online.com Page 1/2 | Product datasheet for ABIN6391940 | 09/12/2023 | Copyright antibodies-online. All rights reserved.

Target Details

. e. get 2 e tee	
Gene:	NPC2
Alternative Name:	NPC2 (NPC2 Products)
Application Details	
Application Notes:	Optimal working dilution should be determined by the investigator.
Restrictions:	For Research Use only
Handling	
Formati	Viral Particlas

Format:	Viral Particles
Storage:	-80 °C
Expiry Date:	12 months
Publications	

Product cited in:

Johnson, Drugan, Miller, Evans: "38" in: , Vol. 1363, Issue Nucleic acids research, pp. 28-39, (1991)